

# PROVOX<sup>®</sup>

## GuideWire

**MD**

**Rx  
ONLY**

**CE** 0413

**STERILE | EO**

42°C  
108°F  
RT  
2°C  
36°F

Fig. 1





Fig. 2

## Prescription information

**CAUTION:** United States Federal law restricts this device to sale, distribution and use by or on order of a physician or a licensed practitioner. The availability of this product without prescription outside the United States may vary from country to country.

## Disclaimer

Atos Medical offers no warranty - neither expressed nor implied - to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

## Patents and trademarks

Provox® is a registered trademark owned by Atos Medical AB, Sweden.

For information about protective rights (e.g. patents), please refer to our web page [www.atosmedical.com/patents](http://www.atosmedical.com/patents).

Compatible Products; Kompatible Produkte; Compatibele producten; Produits compatibles ; Prodotti compatibili; Productos compatibles; Produtos compatíveis; Kompatibla produkter; Kompatible produkter; Kompatible produkter; Yhteensopivat tuotteet; Samhæfar vörur; Ühilduvad tooted; Suderinami gaminiai; Kompatibilni produkty; Kompatibilis termék; Kompatibilné produkty; Združljivi izdelki; Zgodne produkty; Produse compatibile; Kompatibilni proizvodi; Kompatibilni proizvodi; Συμβατά προϊόντα; Совместимые изделия; Uyumlu Ürünler; Համատեղելի արտադրանքներ; Produk yang Serasi; 互換性のある製品; 호환 가능한 제품; 相容產品; מוצרים תואמים; المنتجات المتوافقة

Provox Vega
Provox Vega XtraSeal
Provox2 Voice Prosthesis

## CONTENTS

EN - ENGLISH.....	6
DE - DEUTSCH.....	9
NL - NEDERLANDS.....	13
FR - FRANÇAIS.....	17
IT - ITALIANO.....	21
ES - ESPAÑOL.....	25
PT - PORTUGUÊS.....	29
SV - SVENSKA.....	33
DA - DANSK.....	36
NO - NORSK.....	39
FI - SUOMI.....	43
IS - ÍSLENSKA.....	46
ET - EESTI.....	50
LT - LIETUVIŲ KALBA.....	53
CS - ČESKY.....	57
HU - MAGYAR.....	61
SK - SLOVENČINA.....	64
SL - SLOVENŠČINA.....	68
PL - POLSKI.....	72
RO - ROMÂNĂ.....	76
HR - HRVATSKI.....	79
SR - SRPSKI.....	83
EL - ΕΛΛΗΝΙΚΑ.....	87
TR - TÜRKÇE.....	91
HY - ՀԱՅԵՐԵՆ.....	95
RU - РУССКИЙ.....	99
MS - BAHASA MELAYU.....	103
JA - 日本語.....	107
KO - 한국어.....	111
ZHTW - 繁體中文.....	114
HE - עברית.....	121
AR - العربية.....	125

# 1. Descriptive Information

## 1.1 Indications for use

**Provox GuideWire** is a sterile single use insertion device intended for placement of a sterile Provox indwelling voice prosthesis after total laryngectomy (primary or secondary puncture), or for retrograde replacement of a sterile Provox indwelling voice prosthesis.

## 1.2 Contraindications

Do not use if the patient has anatomical abnormalities, e. g. significant pharyngeal stenosis above the puncture site or severe trismus. Significant pharyngeal stenosis may preclude insertion of the voice prosthesis. Severe trismus may preclude proper protection of the pharyngeal wall during secondary puncture leading to harm of the esophageal tissue.

## 1.3 Description of the device

Provox GuideWire is a sterile, single use device for insertion and replacement of sterile Provox indwelling voice prosthesis.

Provox GuideWire consists of a tube made of a PVC plastic material (Fig. 1.1) and a connector for attachment of a Provox voice prosthesis safety strap made of nylon plastic material (Polyamide) (Fig. 1.2).

Provox GuideWire package contains the following items:

- 1 Provox GuideWire, sterile (Fig. 1),
- 1 Instructions for use Provox GuideWire, non-sterile.

## 1.4 WARNINGS

### Pre-surgery

- **DO NOT** use the product if the package or product is damaged or opened. Unsterile product may cause infection.
- **DO NOT REUSE** and **DO NOT RESTERILIZE** by any method. Intended for single use only. Reuse may cause cross contamination. Cleaning and resterilization may cause structural damage.

### During insertion

- **DO** always proceed slowly and without using excessive force when pulling Provox GuideWire through, failure to follow the warning can lead to tissue damage or subsequent ingestion of parts of Provox GuideWire.
- **DO NOT** proceed if the voice prosthesis gets stuck in pharynx or esophagus, this may lead to tissue damage.

For further information about these events, see section Adverse events.

## 1.5 PRECAUTIONS

Always assess the suitability before using Provox GuideWire in the patient. In cases with lacking suitability, do not use the product.

- **DO** carefully assess patients with bleeding disorders or patients undergoing anticoagulant treatment for the risk of bleeding or hemorrhage prior to placement or replacement of the voice prosthesis.
- **DO** ensure that in case of known or suspected pharyngeal stenosis, dilatation of the pharynx should be performed before using Provox GuideWire.
- **DO** always use aseptic technique during the procedure in order to reduce the risk of infection.
- **DO** ensure that the safety strap of the voice prosthesis is adequately attached to Provox GuideWire.

If proper attachment is not achieved, the voice prosthesis may come loose from Provox GuideWire, causing failure to complete the procedure.

- **DO NOT** load or reload the voice prosthesis if the safety strap has been cut off or is broken.

## 1.6 Adverse events and troubleshooting information

**Aspiration of parts of Provox GuideWire or parts of Provox System** – Accidental aspiration of parts of Provox GuideWire or parts of Provox System may occur. As with any other foreign body, complications from aspiration of a component may cause obstruction or infection. Immediate symptoms may include coughing, wheezing or other abnormal breathing sounds, dyspnea, and respiratory arrest, partial or inadequate air exchange and/or asymmetrical chest movement with respiration. Complications may include pneumonia, atelectasis, bronchitis, lung abscess, bronchopulmonary fistula and asthma. If the patient can breathe, coughing may remove the foreign body. Partial airway obstruction or complete airway obstruction requires immediate intervention for removal of the object. If aspiration of the object is suspected, the object should be located and retrieved endoscopically using a non-toothed grasping forceps.

**Ingestion of parts of Provox GuideWire or parts of Provox System** – Accidental ingestion of parts of Provox GuideWire or parts of Provox System may occur. As with any other foreign body, the symptoms caused by ingestion of parts of Provox GuideWire or parts of Provox System, depend largely on size, location, degree of obstruction (if any) and the length of time it has been present. Ingested components that have remained in the esophagus may be removed by esophagoscopy or observed for a short period of time. The object may pass spontaneously into the stomach. Foreign bodies that pass into the stomach usually pass through the intestinal tract. Surgical removal of foreign bodies from the intestinal tract must be considered when bowel obstruction occurs, bleeding is present, perforation occurs, or the object fails to pass through the intestinal tract.

Spontaneous passage of the object may be awaited for 4-6 days. The patient should be instructed to observe the stools for the ingested object. If the object does not pass spontaneously, or if there are signs of obstruction (fever, vomiting, abdominal pain) a gastroenterologist should be consulted. The device may be retrieved by using a non-toothed grasping forceps.

**Provox GuideWire gets stuck in the pharyngeal mucosal wall** – Provox GuideWire may get stuck in or be interfered by the pharyngeal mucosal wall. With slight pressure, Provox GuideWire will generally bend near the tip and slide upwards towards the pharynx. Stop the procedure if the general bending does not help.

**Hemorrhage/Bleeding of the puncture** – Slight bleeding from the edges of the tracheoesophageal (TE) puncture may occur during intended use of Provox GuideWire and generally resolves spontaneously. Patients on anti-coagulant therapy, however, should be carefully evaluated for the risk of hemorrhage.

**Rupture of the tracheoesophageal tissue** – In case of a rupture of the TE tissue the TE puncture procedure should be abandoned, and the rupture should be sutured immediately. TE puncture should only be repeated after proper healing of the tissues.

**Provox GuideWire gets stuck inside the esophagus** – If Provox GuideWire gets stuck inside the esophagus, do not use excessive force. Attempts to pull out a stuck Provox GuideWire may cause tissue damage and/or breaking Provox GuideWire. The device has to be retrieved with an endoscope.

## 2. Instructions for use

### 2.1 Preparation

Check the integrity of the sterile package. Do not use the product if the package is damaged or opened.

### 2.2 Operating instruction

#### 2.2.1 Primary puncture placement:

1. Create the TE puncture as per your preferred method.
2. Introduce Provox GuideWire in anteroposterior direction through the TE puncture.
3. Guide the tip of Provox GuideWire to the superior margin of the surgical pharyngeal defect.

**CAUTION:** Do not attach the voice prosthesis with a pharynx protector still in situ as the voice prosthesis may get stuck in the protector.

4. Attach the safety strap of the voice prosthesis securely to the tip of Provox GuideWire. Careful traction on Provox GuideWire will bring the strap of the voice prosthesis into the TE puncture.
5. Rotate the voice prosthesis carefully into position with two non-toothed hemostats/forceps holding the tracheal flange.

**CAUTION:** Do not pull the voice prosthesis in place with the GuideWire.

6. When the voice prosthesis is correctly positioned, cut the safety strap off while it remains attached to Provox GuideWire.
  7. Rotate the voice prosthesis into its proper position.
- The laryngectomy procedure can now be completed.

#### 2.2.2 Secondary puncture placement:

1. Create the TE puncture as per your preferred method.
2. Introduce Provox GuideWire in anteroposterior direction through the TE puncture.
3. Guide the tip of Provox Guidewire superiorly into the pharynx and oral cavity and out of the mouth.

**CAUTION:** Do not attach the voice prosthesis with an endoscope still in place as the voice prosthesis may get stuck in the scope.

4. Attach the safety strap of the voice prosthesis securely to the tip of Provox GuideWire. Careful traction on Provox GuideWire will bring the strap of the voice prosthesis into the TE puncture.
5. Rotate the voice prosthesis carefully into position with two non-toothed hemostats/forceps holding the tracheal flange.

**CAUTION:** Do not pull the voice prosthesis in place with Provox GuideWire.

6. When the voice prosthesis is correctly positioned, cut the safety strap off while it remains attached to Provox GuideWire.
7. Rotate the voice prosthesis into its proper position.

#### 2.2.3 Replacement procedure (Fig. 2):

##### **Anterograde vs. retrograde insertion**

If the indwelling Voice Prosthesis cannot be inserted in an anterograde manner, retrograde insertion by use of Provox GuideWire may be performed.

Replacement may be carried out as an outpatient procedure using local anaesthesia. Occasionally, the procedure may be carried out under general anaesthesia upon determination by the physician.

##### **Anterograde removal of the old prosthesis:**

1. Grasp the tracheal flange with a non-toothed hemostat and pull the prosthesis out from the TE puncture.
2. Insert Provox GuideWire through the TE puncture and push Provox GuideWire upwards through the pharynx and out of the mouth.

## **Provox GuideWire insertion through the old voice prosthesis with removal of old voice prosthesis over Provox GuideWire:**

1. While the old voice prosthesis is still in situ, introduce Provox GuideWire through the old prosthesis and direct it through the pharynx and out of the mouth.
2. Grasp the tracheal flange of the old voice prosthesis with a non-toothed hemostat and pull the old voice prosthesis out antegrade and slide it off Provox GuideWire.

### **Placement of indwelling Provox voice prosthesis:**

Provox GuideWire connector is now outside the mouth (Fig. 2.1).

3. Attach and secure the new voice prosthesis to Provox GuideWire Connector. Pull Provox GuideWire in caudal direction. Careful traction on Provox GuideWire brings the strap of the voice prosthesis into the TE-puncture (Fig. 2.2).

**CAUTION:** Do not pull the voice prosthesis in place with the GuideWire.

4. Introduce the tracheal flange of the voice prosthesis into the puncture with the help of a non-toothed hemostat. Then turn the voice prosthesis so that the safety strap is pointing downwards, into the trachea. Let the patient drink and swallow while observing the voice prosthesis, to ensure no leakage occurs. Cut off the safety strap while it is still attached to Provox GuideWire and discard Provox GuideWire (Fig. 2.3).

## **2.3 Cleaning and sterilization**

Provox GuideWire is provided sterile (EO) and intended for single use only and can NOT be cleaned or resterilized.

## **2.4 Disposal**

Always follow medical practice and national requirements regarding biohazard when disposing of a used medical device.

## **3. Reporting**

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

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© Atos Medical AB, Sweden.